

RENA-DRYCON

Haemodialysis Concentrate Mixing System



Just press the mixer button, pour the powder and leave the rest to the Timer & level sensors which fill the required amount of water, mixes and filters a batch of 80 / 160 lit. solution within 10 minutes

Technical data :

Model-	Rena-Drycon
Capacity -	80/ 160 liters
Processing time-	10 mins.
Concentrates-	Acid and Bicarb
Electrical-	230V ac, 15amp, 50Hz
Dimension-	52"height x 26" width

MARKETED BY: SPIN Medlab Private Limited

Office: HIG 1028, 6th Mainroad, Mogappair ERI Scheme, West Mogappair, Chennai - 600037.
+91 44 43235222 | info@spinmedlab.com | www.spinmedlab.com

Manufactured and Marketed by:



ATLANTIC PHARMACEUTICALS

101, Business Point, 349, Western Express Highway, Andheri (East), Mumbai - 400 069.

● Tel: 022-61594343 / 26824313 / 26820950/52. ● Fax : 022-26821627 / 61594311

Email : atlanticpharma@vsnl.com

Website : www.atlanticbiomedical.in