

# RENASTERIL<sup>®</sup>

(Disinfectant for Haemodialysis Machine)

**HOT DISINFECTANT**

Now  
Manufactured in India  
for the First Time  
**FDA APPROVED**

## Application

For use in Haemodialysis Machine with proportional mixing system  
Effective at above 60°C

## Action

Disinfecting and effectively cleaning bio-film, calcium, magnesium deposits to reduce Haemodialysis Machine breakdown

## Composition

21% Citric acid, Malic acid, Lactic acid as adjuvants

## Features

Used for chemo-thermal disinfection of Haemodialysis machine with re-circulation  
Biodegradable

## Package

5-Litre canister

## Storage

Keep container tightly sealed in an upright position at room temperature



MANUFACTURED BY :

**ATLANTIC BIOMEDICAL P. Ltd.**

MARKETED BY : SPIN Medlab Private Limited

Office: HIG 1028, 6th Mainroad, Mogappair ERI Scheme, West Mogappair, Chennai - 600037.  
+91 44 43235222 | info@spinmedlab.com | www.spinmedlab.com

101, Business Point, 349, Western Express Highway, Andheri (East), Mumbai - 400 069.

Tel: 022-61594343 / 26824313 / 26820950/52. Fax : 022-26821627 / 61594311

Email : atlanticpharma@vsnl.com

Website : www.atlanticbiomedical.in