

Laboratory Incubator



Capacity	160 Litre
Chamber Material	Stainless Steel
Temperature Range	+5 To +55 Degree C
Dimensions Exterior	(W 680 x D 691 x H 1008) mm
Dimensions Interior	(W 540xD 510 x H 724) mm
Max. Load Per Shelf	9 Kg
Number of Shelves	4

SPIN Medlab Private Limited

HIG 1028, 6th Main Road, Mogappair ERI Scheme, West Mogappair, Chennai - 600 037.
 +91-9940495550 | +91 44 - 43235222 | info@spinmedlab.com | www.spinmedlab.com
 CIN NO. : U33100TN2019PTC128784.